

# Know and Maintain Your Power Wheelchair



## Maintenance Schedule

### DAILY

#### Frame

Inspect Shroud

#### Electrical System

Charge batteries; Inspect controls, indicators, horn, motor sound, motor disengage lever

### WEEKLY

#### Wheels

Inspect tires

#### Cushion

Inspect cover

### MONTHLY

#### Frame

Inspect welds and suspension

#### Cushion

Inspect cover

#### Wheels

Inspect drive wheels, caster wheels, anti-tip casters; Clean caster axle

#### Supports

Inspect positioning belts, arm, foot, leg, back and other supports

#### Electrical System

Inspect controller box, battery, wiring and electrical connections; Clean power seat functions